



PLASMA EXPRESSOR



Quality, Management, Environment, Health & Safety Standards



CDSCO Certified



USFDA Listed



European CE



ISO 13485:
NABCB/IAF



ISO 9001: IAF



7 Electric Safety
Certifications

INTRODUCTION

The MarkEn Plasma Expressor is expertly designed for the efficient separation of plasma and platelets from blood collection bags. Engineered for reliability, it features a high-tension spring mechanism for quick and seamless plasma transfer, ensuring smooth operation with every use.

KEY FEATURES



Universal Compatibility – Manually operated and suitable for all blood bag types.




Durable & Hygienic Build – Constructed from non-corrosive SS304 stainless steel for longevity and easy cleaning.



Clear Visibility & Precision – Transparent, durable acrylic compression plate for accurate monitoring.



Smooth & Efficient Transfer – Spring-loaded front panel applies uniform pressure, ensuring consistent fluid separation.



TECHNICAL SPECIFICATION	
Make	MarkEn
Material of Body	Corrosion resistance CRCA Steel
Front Panel	Acrylic
Dimension (WxDxH) in mm	675x485x585 (CRCA/GI)
Weight (Net)	3 Kg

CERTIFICATIONS

Quality, Management, Environment, Health & Safety Standards



**50 Products
CDSCO Certified**



USFDA Listed



**8 products tested by a
WHO-accredited laboratory**



European CE



NABCB: ISO 13485



ISO 9001: IAF



ISO 14001: IAF



ISO 45001: IAF



**7 Electric Safety
Certifications**



WHO GMP



**BIFMA Certified,
BIFMA Level 3**



**ASTM- Corrosion
Free Products**



ISO 2409

**EN 13523-26,
ISO 13849-1,
IP 69, IP 65**

& many other...

ALL-IN-ONE SOLUTION FOR

Cold Chain Solution
Medical Equipment
Blood Bank Equipment

Hospital & Scientific Equipment
Hospital Furniture & Infrastructure
Commercial Refrigeration



25+

Years of Excellence



350+

Product Range



500+

Service center



1000+

Committed Engineers



1,00,000+

Satisfied Users

OUR ESTEEMED CUSTOMERS



& MANY MORE...



-  **Regd office & Factory:** 93/1, MIDC Satpur, Nashik- 422007 (IND)
-  **Corporate office:** 410, Maker Chambers V, Nariman Point, Mumbai - 400021
-  **New Delhi office** - DLF Tower, Jasola | **Chandigarh Office** - CCC, VIP Road
-  **Customer Care No.**+91 91 3030 5959
-  **www.markenworld.com**  **info@markenworld.com**



Authorized Distributor Details