



PLASMA EXPRESSOR



Quality, Management, Environment, Health & Safety Standards



CDSCO MD 5 License



USFDA Listed



European CE



ISO 9001: 2015
IAF



ISO 13485:2016
NABCB/IAF



ISO 45001 : 2018
IAF



7 Electric Safety Certifications
IEC 60335-1 and IEC 60335-2-24, IEC 61010-1,
IEC 61010-2-020 and 61010-2-101, IEC 61326-1,
IEC 61326-2-6.

INTRODUCTION

The MarkEn Plasma Expressor is expertly designed for the efficient separation of plasma and platelets from blood collection bags. Engineered for reliability, it features a high-tension spring mechanism for quick and seamless plasma transfer, ensuring smooth operation with every use.

KEY FEATURES



Universal Compatibility – Manually operated and suitable for all blood bag types.




Durable & Hygienic Build – Constructed from non-corrosive SS304 stainless steel for longevity and easy cleaning.



Clear Visibility & Precision – Transparent, durable acrylic compression plate for accurate monitoring.



Smooth & Efficient Transfer – Spring-loaded front panel applies uniform pressure, ensuring consistent fluid separation.



TECHNICAL SPECIFICATION

Make	MarkEn
Material of Body	Corrosion resistance CRCA Steel
Front Panel	Acrylic
Dimension (WxDxH) in mm	675x485x585 (CRCA/GI)
Weight (Net)	3 Kg

ALL-IN-ONE SOLUTION FOR

Cold Chain Solution
Medical Equipment
Blood Bank Equipment

Hospital & Scientific Equipment
Hospital Furniture & Infrastructure
Commercial Refrigeration



25+

Years of Excellence



400+

Product Range



500+

Service center



1000+

Committed Engineers



1,00,000+

Satisfied Users

CERTIFICATION



10 products tested by a WHO-PQS accredited laboratory



ISO 14001: IAF



WHO GMP



BIFMA Certified, BIFMA Level 3



ASTM- Corrosion Free Products



ISO 2409

EN 13523-26,
ISO 13849-1,
IP 69, IP 65

& MANY OTHER...

OUR ESTEEMED CUSTOMERS



& MANY MORE...



Regd office & Factory: 93/1, MIDC Satpur, Nashik- 422007 (IND)

New Delhi office - 05, A-44, VDS Tower, Sector 2, Noida- 201301. UP. India

Customer Care No.+91 91 3030 5959

www.markenworld.com info@markenworld.com

